

PROJECT

DATE

TITLE

Continued from page

4) streak out strains : DH10B(DE3)
 DH10B pLysS
 BL21 AI
 BL21(DE3) pLysS
 DH10B? / TOP10?

5) prepare TUCOI template

6) order primers

20180619 Cloning successful. - pAAG169 is PET9a GG-GFP vector
 pAAG170 is PET9a-6xH-TEV-LacII-331

Made frozen stock of SAAG170a (DH10B).

Transformed pAAG170 →
 1) BL21 AI
 2) TOP10
 3) DH10B pLysS
 4) DH10B(DE3)
 5) BL21(DE3) pLysS

Began O/N cultures of 20 ml each of SAAG114b, SAAG166b for expression.
 Began 2 ml O/N cultures of all above empty strains as well.

20180620 Made frozen stocks of empty e. coli strains.
 Transformation plates of pAAG170 in all strains ⇒
 lots of colonies.
 ⇒ ready for expression - stored plates at 4°C.

S9610 MBPα Y197A, ARm3 protein expression cultures : 1L each.

- Subcultured 20 ml → 1L LB/Kan/Cm O/N culture
 + 5 ml 40% glucose for 114b
- Grew 2.5 hrs at 37°C / 225 rpm
- Induced with 500 μM IPTG and reduced temperature to 16°C.
- Grew cultures O/N.

Plan for crystallography : - dialyze into pH 8.3-8.5 buffer to set up
 Mosquito trays / 15-well S9610 trays.
 - optimize PEG3350 / Na Acetate conditions
 from previous Mosquito hits (after dialysis
 into pH 6.0 MES.)
 ⇒ JCSG1 / CB

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE

Continued from page

20180621

Spun down SAAG114b, SAAG146b expression cultures at 6000 x g / 4°C / 20 min. Froze pellets.


Thawed 114b pellet on ice (BL21(DE3) pLysS, pET28-MBP-Y197A). Added 2.25 mL DNase I, 225 mL MgCl₂, 45 mL MnCl₂, 4.5 mL CaCl₂, and protease inhibitor to 40 ml 50 mM Tris, 150 mM NaCl, pH 8.0 buffer. Resuspended pellet in lysis buffer by vortexing for 15-20 min.

Lysed cells by 2x microfluidizer. Spun down insoluble fraction at 27000 x g, 30 min, 4°C. Decanted supernatant and stored on 4°C for 1 hour. Spun down at 3700 x g / 4°C / 10 min. and decanted again. Passed through freshly regenerated amylose column. Eluted in 14 ml Tris buffer + 30 mM maltose. Concentrated down to 1 mL in 10 kDa MWCO spin column, then applied to S200 for SEC overnight.

20180622

MBPY197A SEC: 3 peaks

equilibrate_elute_equilibrate_SEC_Y197A_20180621 001



Peak Table - UV 1280

Peak	Retention ml	Area mAU	Area %	Ext coeff. mg ml ⁻¹ cm ⁻¹	Fraction(s)	Volume ml	Conductivity mS/cm
Peak A	69.330	322.5	15.5		4.B.8 - 4.C.3	16.023	15.41
Peak B	75.960	1238	59.51		4.C.5 - 4.C.8	7.999	15.41
Peak C	86.037	519.8	24.99		4.C.10 - 4.D.3	11.998	15.42

Thawed ARm3 pellet on ice and resolubilized in lysis buffer. Lysed by microfluidizer and spun down the insoluble fraction.

Continued to page

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE

Continued from page

20180623 ran gel on purified proteins. Destaining o/n.
 Prepared dialysis buffer: 25 mM Tris, 150 mM NaCl, pH 8.0.
 Began o/n cultures: (1) BL21 AI / pAAG170
 (2) BL21(DE3) pLys S / pAAG170
 (3) DH10B(DE3) pLys S / pAAG170 } 21 ml ea.

20180624 Made frozen stocks:
 (1) SAAG170b : BL21 AI
 (2) SAAG170c : BL(DE3) pLys S
 (3) SAAG170d : DH10B(DE3) pLys S

Began protein expression cultures:

- subcultured o/n 20 ml → 1L LB + cm/kan for (2,3)
 + Kan for (1)
 + 5 ml 40% glucose
- grew at 37°C / 225 rpm for 2-2.5 hours until OD₆₀₀ ~ 0.4
- induced with 10 ml 20% arabinose (final conc. 0.2%) for (1)
 and with 1mM IPTG for (2,3)
- reduced temperature to 16°C at induction and grew o/n.

OD₆₀₀ table —

	1 hour	1.5	2	2.5
(1)	0.13	0.29	0.49	
(2)	0.12	0.24	0.42	
(3)	0.08	0.13	0.20	0.45

SDS-PAGE

- 1 load MBPY197A SEC
- 2 SEC Y197A peak 1
- 3 SEC Y197A peak 2
- 4 SEC Y197A peak 3
- 5 load ARM3 HTQ
- 6 SEC ARM3 peak
- 7 old sample Y197A
- 8 old sample ARM3.

Combined fractions and dialyzed
 at 4°C, spinning, o/n.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

PROJECT

DATE

TITLE

Continued from page

20180625

Spun down o/N cultures, 6000 x g / 4°C / 20 min. Freeze pellets at -80°C.
 Took proteins ARm3, 59 G10 MBP & Y197A out of dialysis and concentrated down to ~200 µl volumes.

Stocks: ARm3 : 2.2 mM / 38.9 mg/ml
 Y197A : 970 µM / 39.5 mg/ml

Prepared protein solution for Mosquito trays (pH 8.0): 200 µl total
 15.45 µl ARm3
 35.0 µl Y197A
 39.0 µl FPP (from 10 mM stock)
 115.55 µl buffer (Tris / NaCl pH 8.0)

Buffer-exchanged into 25 mM MES, 200 mM NaCl, pH 6.0 for JCSG1 / C8 optimization trays. Conditions optimized around: 0.2 M Na Acet
 20% w/v PEG3350
 ↳ set up one 15-well tray w/ 3 concentrations of protein / FPP

Also set up 6 x 15-well trays with 3 concentrations in each well of pH 8.0 ARm3 / Y197A in HIGH PEG 59 G10 conditions — 3x no maltose, and 3x 170 µM maltose.

Protein final concentrations: 170 µM, 85 µM, 70 µM.
 FPP at 10x protein concentration in each solution.

Stored all trays in crystallography room.

20180627 OCTET experiments: large concentration range
 ARm3 / avi-Y197A BT
 + 200 µM FPP.

Bridge stocks: - non-BT avi-Y197A @ 0.9 mg/ml ~ 21 µM (large epit tube)
 - Y197A-BT @ 1.9 mg/ml ~ 40 µM (2ml tube)
 ↳ diluted 12 µl → 238.6 µl buffer for 200 nM experimental conc. (200 µl / well)

↳ buffer: 50 µM Tris, 150 mM NaCl, pH 8.0
 - ARm3, 2.2 mM fresh conc. stock solution.

↳ diluted 5.5 µl (2.2 mM) → 49.5 µl buffer
 for 220 µM ARm3 stock
 ↳ diluted 1.1 µl (220 µM) → 9.9 µl buffer for 22 µM ARm3 stock.

PROPRIETARY INFORMATION

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

TITLE

Continued from page

experiment: Need 55 μ L 220 μ M ARM3, 11 mL 22 μ M ARM3, 2.4 mL 220 μ M avi-Y197A-BT, 36 μ L 10 mM FPP

	sample conc.	sample ID	sample volume	buffer volume
A1	0 μ M	22 μ M ARM3	0 mL	196.0 μ L
B1	0.1 μ M		0.9 μ L	195.1 μ L
C1	0.2 μ M		1.8 μ L	194.2 μ L
D1	0.3 μ M		2.7 μ L	193.3 μ L
E1	0.5 μ M		4.5 μ L	191.5 μ L
F1	1.0 μ M	220 μ M ARM3	0.9 μ L	195.1 μ L
G1	1.5 μ M		1.4 μ L	194.6 μ L
H1	3.0 μ M		2.7 μ L	193.3 μ L
A2	6.0 μ M		5.5 μ L	190.5 μ L
B2	8.0 μ M		7.3 μ L	188.7 μ L
C2	12.0 μ M		10.9 μ L	185.1 μ L
D2	24.0 μ M		21.8 μ L	174.2 μ L

★ +4 μ L FPP/well from 10 mM stock

PROTOCOL:

1. RINGE / 60 s
2. LOAD / 300 s
3. RINSE / 60 s
4. BASELINE / 60 s
5. ASSOCIATION / 120 s
6. DISSOCIATION / 240 s

Subtracted 0 μ M AR well from each sample.

Local / partial fitting
Req. fit for steady-state analysis.

Ran the experiment 4x.

1:1 MODEL ~ the fit wasn't so great for any of the trials.

Continued to page




DATE

SIGNATURE


DATE

PROPRIETARY INFORMATION


DISCLOSED TO AND UNDERSTOOD BY

Sensor Type: SA (Streptavidin) - by Conc. (μM)Sensor Type: SA (Streptavidin) (SS)
KD=1.5E-07 Rmax=0.2787 R²=0.7648Sensor Type: SA (Streptavidin) - by Conc. (μM)

Sensor Type: SA (Streptavidin) (Iso-Affinity)



Sensor Type: SA (Streptavidin) (X-Y)



#1 - ALL DATA.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION